


Module Treatment Simulation 2D

Module Treatment Simulation 2D can be used to simulate orthodontic and/or surgical treatment measures on a lateral cephalometric image in combination with an associated lateral photo in their effect on cephalometric values and soft tissue profile changes.

First Steps

1. Open module wit lateral ceph xray and (optionally) a corresponding lateral photo
2. Define segments for simulation
3. Select method and parameters for soft tissue simulation
4. Displace / rotate segments
5. Analyze cephalometric changes (graphically, numerically)
6. Perform soft tissue simulation: compare before and after
7. Save planning as a project
8. Save modified images (lateral ceph, lateral photo) as findings

Module user interface



Links



From:
<http://onyxwiki.net/> - [OnyxCeph^{3™} Wiki]



Permanent link:
http://onyxwiki.net/doku.php?id=en:treatment_simulation_2d&rev=1584542956

Last update: **2020/03/18 15:49**