


Panel Objektliste

Zu jeder 3D-Szene (z.B. im 3D-Modul, 3D-Vollbild) kann eine Objektliste angezeigt werden. In ihr sind alle Teilobjekte des aktuell geladenen Datensatzes aufgelistet. Einzelne oder mehrere Listeneinträge können mit linkem bzw. mit [STRG] + linkem Mausklick ausgewählt werden.

Ist kein Teilobjekt ausgewählt, beziehen sich Objektmanipulationen sowie Änderungen an den [Objekt-](#) und [Ansichtsoptionen](#) auf alle Teilobjekte der Liste.

Sind ein oder mehrere Teilobjekte ausgewählt, beziehen sich Objektmanipulationen sowie Änderungen an den Objekt- und Ansichtsoptionen auf alle ausgewählten Teilobjekte der Liste.




Kontextmenü:

Im Kontextmenü der Objektliste sind verschiedene Einträge zu Vorauswahl, Sichtbarkeit, Netztyp, Einfärbung u.a. vorhanden. Die im Kontextmenü aufgeführten Einträge können sich in einzelnen 3D-Modulen je nach Modulcharakteristik unterscheiden.

Eintrag	Bedeutung

Statuszeile:

Bei Doppelklick auf die Statuszeile am unteren Rand der Objektliste werden die Ergebnisse einer topologischen Analyse des Oberflächennetzes für das in der Liste ausgewählten Teilobjektes angezeigt.



The screenshot displays the OnyxCeph3 software interface. On the left, a 3D model of a dental arch is shown with a blue retainer band. A central window titled "OnyxCeph^{3™}" provides technical data:

- Vertices: 2096, Faces: 3140
- Possible pathologies:
 - Isolated vertices: 0
 - Non-manifold vertices: 0
 - Isolated Faces: 0
 - Faces topo. collapsed to point: 0
 - Faces topo. collapsed to edge: 0
 - Faces geom. collapsed to point: 0
 - Faces geom. collapsed to edge: 0
 - Faces without area: 0
 - Non-manifold edges: 0
 - Edges with wrong orientation: 0
- Boundaries:
 - Number of boundaries: 4 (262/262/262/262 Vertices)
 - Number of boundary errors: 0
- Shells:
 - Number of shells: 3 (262/262/1572 Vertices)
- Quality:
 - Triangles with minimum angle > 32°: 1640
 - Triangles with minimum angle > 16°: 1405
 - Triangles with minimum angle > 8°: 90
 - Triangles with minimum angle > 4°: 5
 - Triangles with minimum angle > 2°: 0
 - Triangles with minimum angle > 1°: 0
 - Triangles with minimum angle < 1°: 0
 - Invalid triangles: 0
 - (0% of all faces are of bad quality (< 4°) or invalid.)
- OC:T=103,G=0

An "Objektliste" panel on the right lists objects with checkboxes. The "Retainer-Test Maxilla" and "Retainer-Test Mandibula" are checked. The bottom of the panel shows "1 Objekt (•:2.096 Δ:3.140)". A red arrow points to the "Retainer-Test Maxilla" entry with the text "2 x Left".

From:
<http://onyxwiki.net/> - [OnyxCeph^{3™} Wiki]

Permanent link:
http://onyxwiki.net/doku.php?id=panel_objectlist&rev=1446329382

Last update: 2015/10/31 23:09

