


Hand over Master Data



The master data of the currently selected patient are written (optionally according to UTF-8 coding) to the defined text file or to the Windows clipboard.

If the call is made from the Patient tab, the patient ID, surname, first name, birthday and gender are transferred.

```
Onyx Data
10014
Doe
John
19800101
M
```

If the Images tab is called, the subdirectory of the image file in server folder \\Pictures, the image code (Synccode) and the image acquisition date of the focused image data set assigned to the image in the session overview are output in addition to the master data.

```
Onyx Data
10014
Doe
John
19800101
M
10014\100141AFU.ii3
100141AFU
20210126
```

From:
<http://onyxwiki.net/> - [OnyxCeph^{3™} Wiki]



Permanent link:
<http://onyxwiki.net/doku.php?id=en:homasterdata>

Last update: **2021/08/13 16:13**