

# Documentation

<a href="#">Features and Prices</a>
<a href="#">Setup Instructions</a>
<a href="#">Version Overview</a>
<a href="#">License Types</a>
<a href="#">Onyx Communication Options</a>
<a href="#">Move Database Server 2.7 Instruction</a>
<a href="#">Move Database Server 3.x Instruction</a>
<a href="#">Requirements for Customized Analysis</a>
<a href="#">Instructions For Use</a>
<a href="#">Spezification Face Scanner</a>
<a href="#">Spezification Model Scanner S600</a>
<a href="#">Spezification Model Scanner orthoX</a>
<a href="#">Virtual Treatment Planning</a>
<a href="#">Virtual Orthodontic Treatment Planning</a>
<a href="#">Digital i/d Bonding</a>
<a href="#">3D Scanning - and then what?</a>

## Requirements for Customized Analysis

The following information is required to create a customized analysis:

0. Tracing type of the analysis (Lateral Ceph, Pano, Profile Photo, Frontal Foto, ...)
1. Name of the analysis
2. List of all numeric results to be calculated, incl.
  - 2.1 Short name of the result variable
  - 2.2 Long name of the result variable
  - 2.3 Definition for the result variable incl. no. of digits after decimal point
  - 2.4 Clinical norm (if available)<sup>1</sup>
  - 2.5 Standard deviation for the norm (if available)<sup>1</sup>
  - 2.6 Verbal result description for result interbvals (if available)
3. A sketch for the graphical tracing
4. A sketch for the print form

<sup>1</sup> if needed: norms depending on age and gender; also specify ethnical group, if other than caucasian; It is also accepted to refer to an already implemented analysis – however, make sure that results put together from different existing analyses are harmonized (see OnyxCeph<sup>3™</sup> menu item Help |Findings lexicon).

For creating a customized analysis a flat rate of EUR 290 plus VAT (if applicable) is charged one-off. The analysis is delivered as a patch file to be imported via menu item Options|Install Patch.

From:

<http://onyxwiki.net/> - [OnyxCeph<sup>3</sup>™ Wiki]

Permanent link:

[http://onyxwiki.net/doku.php?id=en:support\\_manuals&rev=1588071509](http://onyxwiki.net/doku.php?id=en:support_manuals&rev=1588071509)

Last update: **2020/04/28 12:58**

