

System Options - Tab |Upload Data / Container|

On system options tab |Upload Data|, all required pre-settings for the communication to web servers via FTP or HTTP protocol are configured to upload or download [Container](#)-, report- or [WebViewer](#) files.

System Options

Patient Data

Image Types

Sessions

Galleries

Printouts

Forms

Display

Interfaces

Third Party Tools

Event

Upload Data

Email

ContainerReportImport3D-Web-Viewer

Host:

Connection Security:

None

User Name:

User Password:

.....

Host Directory:

Passive Mode:

☒

Download - Protocol / Data Address:

https://

onyx.orthorobot.com/container

?

Help

Close

[Click on tab to switch to related topic]

Server Settings Data Upload




The upload has to be configured in [System Options](#) on tab [Upload Data](#) > 3D-Web-Viewer before it becomes available as menu item.

Entry	Description	Example
Host	Your provider for the FTP-webospace	123.456.789.012
Connection Security	Upload protocol	FTPS
User Name	The user name für the FTP access	onyxUser
User Password	The password for the FTP access	••••
Host Directory	The directory of the data relative to the root directory of the user	container
Passive Mode	Passive transfere mode activation	✓

Download Protocol / Data Address

Protocol (FTP | http: | https:) and web address (URL) to download from the server location defined in section 1 above.

Links

	Container-Exchange
	Workflow Container Export
	Workflow Container Import

From:

<http://onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:

http://onyxwiki.net/doku.php?id=en:tab_uploaddata

Last update:

2023/04/14 15:10

