

## #3429 Aligner\_3D: Option Scan Refinement

---

Since release 3.2.7 (77) in module Aligner\_3D, an additional panel is available bottom down in the module window middle panel which has 3 separate buttons. The functionality linked to this panel can be used to adjust a calculated aligner step by the actual treatment situation using an intermediate scan to refine the current step and therefore, the subsequent planning.

**see also:** [Panel \[Scan-Refinement\]](#)

From:

<http://onyxwiki.net/> - [OnyxCeph<sup>3</sup>™ Wiki]

Permanent link:

<http://onyxwiki.net/doku.php?id=en:3429&rev=1445805960>

Last update: **2015/10/25 21:46**

