

Module Bonding Jigs 3D




Program version	BASIC	2D PRO	3D PRO	LAB	OMS
Availability	✗	✗			

Design transfer trays for bracket positions planned in modules such as [FA_Bonding](#) or [Wire_Bonding](#). Specialized adapter objects are required.

First Steps

1. Open Module with bracket-attached Model
2. Select Crowns
3. Assign adapter and Jig parameters to crowns
4. Construct Jigs individually or connected
5. Arrange Jigs for manufacturing / export
6. Save planning as project and/or as finding

Module user interface



Links





From:

<http://onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:

http://onyxwiki.net/doku.php?id=en:bonding_jigs_3d&rev=1586011802

Last update: **2020/04/04 16:50**

