

Module Bonding Jigs 3D




Program version	BASIC	2D PRO	3D PRO	LAB	OMS
Availability	×	×			

In module Bonding Jigs 3D module, occlusal transfer jigs can be constructed to transference the virtually planned bracket positions using metallic standard bracket adapters.

First Steps

- 1. Open Module with bracket-attached Model
- 2. Select Crowns
- 3. Assign adapter and Jig parameters to crowns
- 4. Construct Jigs individually or connected
- 5. Arrange Jigs for manufacturing / export
- 6. Save planning as project and/or as finding

Module user interface



## Links



From:

<http://onyxwiki.net/> - **[OnyxCeph<sup>3</sup>™ Wiki]**

Permanent link:

[http://onyxwiki.net/doku.php?id=en:bonding\\_jigs\\_3d&rev=1587498161](http://onyxwiki.net/doku.php?id=en:bonding_jigs_3d&rev=1587498161)

Last update: **2020/04/21 21:42**

