Client Setup

TP Onyx Admin

Note: This instruction is needed if a new client in the network of an existing clinic server has to be installed

1. Install Client

Client installation on new workstation requires connection to the server via local IP network:

W Client setup

2. Establish Server Connection

Manual configuration of server connection:

W Server connection

3. Create Onyx_PDF Folder

Create folder C:\Onyx_PDF on client.

TP Onyx Admin

From: https://onyxwiki.net/ - [OnyxCeph³™ Wiki]

Permanent link: https://onyxwiki.net/doku.php?id=en:moh:tp_admin5

Last update: 2022/03/15 13:02

