


Tracing Options

Using panel button [Options], various settings on individual tabs with regard to graphics overlays on 2D and 3D image windows can be made.


Tab |Input Window|



[Click on tab opens related wiki page]


Panel Font

Here, font, font color, and background color of reference point labels (point labels, tooltips) can be defined. The settings affect the local client



Panel Lines

Here, for line types Style 1 ... Style 6 selected in the graphic on the left or in the pulldown, size, color and type of the line can be defined.
The settings affect the local client.



From:
<http://onyxwiki.net/> - [OnyxCeph^{3™} Wiki]



Permanent link:
http://onyxwiki.net/doku.php?id=en:to_input

Last update: **2021/10/14 15:35**