

# Durchzeichnungstypen 3D

OnyxCeph<sup>3™</sup> supports the following tracing types 3D:

Type
Head 3D
Cast 3D Deciduous
Cast 3D Mixed
Cast 3D Permanent
Not Defined 3D

Each 3D image type can have assigned multiple tracing types 2D and 3D. Thus, in module [Digitize 3D](#), after selecting the tracing type to be used, all analyses predefined for this tracing type are offered for selection and application. The tracing type assignment can be modified via menu [Options|System Options](#) on tab [Image Types 3D](#).

From:

<http://onyxwiki.net/> - [OnyxCeph<sup>3™</sup> Wiki]

Permanent link:

[http://onyxwiki.net/doku.php?id=en:tracingtypes\\_3d](http://onyxwiki.net/doku.php?id=en:tracingtypes_3d)

Last update: **2020/07/23 11:07**

