


V.T.O.3D

First Steps

1. Select useful wizard steps with right and apply with left click
2. Adjust bow and single teeth to visual treatment objective as needed
3. Show inter tooth and occlusal distances and adjust planning as needed
4. Save intermediate steps and final result als project
5. Save as new finding for Aligner or Wire-Bonding

Module user interface



References

Screenshots	
3D Views	
Videos	
Links	Request Manual

From:
<http://onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:
<http://onyxwiki.net/doku.php?id=en:vto3d&rev=1563013921>

Last update: **2019/07/13 12:32**

