

Launching Database Server

In standard operation mode, Onyx database server program [OnyDBServer.exe](#) should run as Windows™ service. The service can be launched and configured either by the Windows™ control panel or (recommended) in program window [Onyx3Config](#) by button [Start Server as Program].

If Onyx database server program settings [OnyDBServer.exe](#) should be modified, the program has to be launched as application by button [Start Server as Program] and then opened from the task bar system tray icon.



References

Screenshots	
3D Views	
Videos	
Links	OnyDBServer - Onyx3Config

From: <http://onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link: http://onyxwiki.net/doku.php?id=en:workflow_0016

Last update: 2016/02/05 13:11

