


#3835: Module Retainer 3D: Distance Measurement Retainer-Opposite Jaw

Since Release 3.2.61(131) in module [Retainer 3D](#) by the new panel Distance, an option was added to measure the perpendicular distance between a point placed manually onto the opposite jaw and the nearest upper resp. lower retainer plane.



From:
<http://onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:
<http://onyxwiki.net/doku.php?id=en:3385>

Last update: 2016/07/06 12:56

