

Module Edit 2D



Program version	BASIC	2D PRO	3D PRO	LAB	OMS
Availability	✗				

not included

included by default


optionally included

Module Edit 2D is used to edit / manipulate 2D images by appying different image processing algorithms. For further use, the result of the processing must be saved as a copy.

First Steps

1. Mark region
2. Select filter
3. Activate processing
4. Save result as finding

Module user interface



Links

n/a

From:
<https://onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:
https://onyxwiki.net/doku.php?id=en:edit_2d&rev=1676908875

Last update: **2023/02/20 17:01**

