


# Interface Configuration Digital Xray: Orthometric



Parameter	Setting
[Support Activated]	Activates import icon button in module <a href="#">Add Image 2D</a>
Datenbank [...]	Pfad zur Orthometric Datenbank
Bilder [...]	Pfad zum Orthometric Bildpfad
[Test]	Zugriff auf interne DB testen

From:

<http://onyxwiki.net/> - [OnyxCeph<sup>3™</sup> Wiki]

Permanent link:

<http://onyxwiki.net/doku.php?id=en:orthometric>

Last update: **2022/06/27 07:36**

