

# PMS Interface TK2 (Denmark)

Practice management program from Denmark

## Interface for patient transfer TK2 to OnyxCeph<sup>3</sup>™

TK2 supports Onyx command line interface via third party programs to be configured in TK2.

The screenshot shows the 'Programindstillinger' window with the 'Røntgen' tab selected. Under 'Røntgen', the 'OnyxCeph' sub-tab is active. The 'Aktiveret' checkbox is checked. The configuration fields are as follows:

- Navn: OnyxCeph
- Sti til program: C:\Programmer\Image Instruments\Onyx3Client (Find.. button)
- Navn på program: OnyxCeph.exe (Find.. button)
- Sti til data: (Find.. button)
- Sti til logo: \\tand-adm-02\icon\$ (Find.. button)
- Logo: Onyx\_Logo\_28.bmp (Find.. button)

Below these fields is a 'Parameterliste' (Parameter list) with the following items:

- Patient.Patient (selected)
- Patient.Aktiv
- Patient.AktivDato
- Patient.PasivDato
- Patient.PersonNum
- Patient.KaldeNavn

To the right of the parameter list is the 'Opætning af komando linie' (Command line setup) section, which contains the following command:

```
"C:\Programmer\Image Instruments\Onyx3Client\OnyxCeph.exe" -PatId "<Patient.PersonNum>" -Firstname "<Patient.Fornavn>" -Lastname "<Patient.Efternavn>"
```

A warning message below the command line states: 'Bør kun ændres i samråd med DS/ DATA!'.

From:

<http://www.onyxwiki.net/> - [OnyxCeph<sup>3</sup>™ Wiki]

Permanent link:

<http://www.onyxwiki.net/doku.php?id=en:tk2>

Last update: 2022/02/04 10:07

