

Window [System Options]

Tab	Settings / Konfiguration ...
1 Tab Patient Data	Patient--thumbnails, VDDS settings, Markers, Misc ...
2 Tab Image Types	2D and 3D Image Type Settings
3 Tab Sessions	2D View Settings 1/2/3, 3D View Settings
4 Tab Galleries	Galleries, Presentations, Sequences, Media, View Option ...
5 Tab Printouts	Printform Templates, Letter Heads, Printer Preconfiguration ...
6 Tab Forms	DHTML-Forms by Categories, Password Forms in Fullscreen Mode ...
7 Tab Digitize	View Options, User Analyses, Use Drawing Templates, Customization Landmark Sequence ...
8 Tab Interfaces	Configuration External Image Sources, External Documents Archive ...
9 Tab Third Party Tools	Configuration 3rd Party Tools
10 Tab Event	Event Settings definieren
11 Tab Upload Data	FTP-/WEB-Server Upload Settings for Containers, Reports, Webviews
12 Tab E-Mail	Email Client and Attachments konfigurieren

From:
<http://onyxwiki.net/> - [OnyxCeph^{3™} Wiki]

Permanent link:
http://onyxwiki.net/doku.php?id=en:wdw_systemoptions

Last update: **2021/01/14 14:43**

