

Launch the archive process by hidden short cuts CTRL+doubleclick on label "Image Folder" in the OnyxDBServer window:

Archive Patient

Mark patient record for next archiving run or undo archiving this way to recover the archived record. Patient record will stay in the patient list if display of archived patient records is activated

Settings Menu:

- ☒ Show Archived Patients
- ☐ Sort Active Column
- ☒ Narrow Down Search Results
- ☐ Live Request
- ☒ Allow Combined Markers
- ☐ Cache List of Patients
- ☒ Update Cache Automatically
- ☐ Memorize Search Items
- ☒ Fit Columns

| Patient's ID | Last Name | First Name | Gender | Model |
|--------------|-----------|--------------|--------|--------|
| 10000 | Demo | Daniela | Female | MN 550 |
| 10101 | Demo | Daniela | Female | |
| 10228 | Demo | Daniela | Female | Doku |
| 10432 | Demo | Daniela | Female | |
| 10000A | Demo | Daniela 이후 변 | Female | MN 550 |
| 10240 | Demo | Leipzig | Male | |

From:

<https://onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:

https://onyxwiki.net/doku.php?id=img_4243_3

Last update: 2021/09/16 17:28

