

Eintrag GingivaShowMovingPatchesVersion

Syntax

```
[Aligner]  
GingivaShowMovingPatchesVersion=value    |    value = 0,1
```

```
[VirtualSetup]  
GingivaShowMovingPatchesVersion=value    |    value = 0,1
```

Beschreibung

Aktiviert Zahnfleischberechnungsmethode "Fixiert"

From:
<http://onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:
http://onyxwiki.net/doku.php?id=ini_aligner_gingivashowmovingpatchesversion

Last update: **2020/03/03 08:18**

