

Client-Programm OnyxReport



Um die in OnyxCeph³™ vorhandenen Informationen auch für Personen verfügbar zu machen, die selbst über keine Softwarelizenz verfügen, können innerhalb des Hauptprogramms Reports zusammengestellt und als Datei lokal gespeichert oder auf einen FTP- oder Webserver hochgeladen werden.

Programmoberfläche

Report

Patient

FB2354 - Roots, Tooth
4178 - Klüber, Johannes

Daten

☒ Send Anonymously

Options

Recipient's language
English

☐ Password

☒ Send Email with FTP-Upload

Images **Slides**

No.	Type	Date of Record	Image
☒ 1	CT-HeadScan	03.03.2022	☒
☒ 2	Lateral Right	03.11.2021	☒
☒ 3	Cast Deciduous 3D	23.11.2021	☒
☒ 4	Cast Mixed 3D	23.11.2021	☒

Comment about Patient

Add patient specific comments here

Sender

Image Instruments

Message for recipient

Add a message for the recipient here

SAMPERMANS IOS INTERNATIONAL ORTHODONTIC SEMINARS

Help Preview Close

Verweise



OnyxReport verwenden



OnyxReport

From:

<http://onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:

http://onyxwiki.net/doku.php?id=program_onyxreport

Last update: **2023/02/16 22:22**

