

Schnittstellenkonfiguration Scanner: ZirkonZahn S300 / S600

Model- i/o-Scanner

CS_3X00	CADstar	DentalDirekt	Imetric
CEREC Ortho	SHINING 3D	DOF	iTero
smart optics	orthoX	R+K	ZirkonZahn
True Definition			



☒ **ZirkonZahn - Support Activated**

ZirkonZahn Data Folder
Z:\ZZ_Daten

ZirkonZahn Configuration
[Empty dropdown menu]

Rotation

1.	Y-Axis	180,0
2.	X-Axis	0,0
3.	X-Axis	0,0

? Help ✓ OK ✕ Cancel

Parameter	Einstellung
-----------	-------------

From:

<https://www.onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:

<https://www.onyxwiki.net/doku.php?id=zirkonzahn>

Last update: 2023/04/28 17:15



